

Applications



Emergency

Monitoring electrolytes in critically ill patients due to poisoning, coma and convulsions



ICU

Monitoring electrolytes in critically ill patients



Surgery

Monitoring electrolytes during operation



Anesthesiology

Monitoring electrolytes during surgical anesthesia (preoperative, intraoperative and postoperative)



Dermatology

Monitoring electrolytes in patients with severe trauma, burn and scald



Gastroenterology

Monitoring electrolytes in patients with diarrhea and vomiting accompanied by coma



Nephrology

Monitoring electrolytes for patients in dialysis ward



General Physicians' office

Primary Medical Electrolytes test

How to use



1. Scan the barcode on the cartridge pouch and take out the cartridge.



2. Fill the inlet with sample to the fill mark and slide the cap to seal the inlet.



3. Insert the cartridge into the analyzer until it clicks. Wait for the results.

Electrolyte Analyzer



Electrolyte Analyzer

Advanced Features

- Concurrent testing K⁺, Na⁺, Cl⁻, iCa²⁺, and iMg²⁺
- Ionized Magnesium (iMg²⁺) testing.
- Dry electrochemical technology with a patented pellet structure design.
- Exclusive molecular probes and algorithms.

Efficiency and Accuracy

- Swift calibration and testing in under 5 minutes.
- Built-in calibration solution.

User-friendly Experience

- Whole blood testing for instant results.
- Minimum sample volume requirement (100 µL).
- No fluidic pathways for hassle-free operation and minimal maintenance.
- Room temperature stored test cartridges (up to 8 months).

Portability and Data Storage

- Compact and portable design for on-the-go use.
- Over 50 consecutive tests on a single charge with the built-in lithium battery.
- Built-in barcode scanner and thermal printer.
- Stores over 99,999 records.
- Connection to LIS/HIS.



Performance

Parameter	Accuracy	Precision	Stability	Test Range (mmol/L)
K ⁺	±3.0%	≤1.5%	≤2.0%	1.0-15.0
Na ⁺	±3.0%	≤1.5%	≤2.0%	100-200
Cl ⁻	±3.0%	≤1.5%	≤2.0%	65-160
iCa ²⁺	±5.0% ±0.05% mmol/L	≤1.5%	≤3.0%	0.25-4.00
iMg ²⁺	±5.0% ±0.05% mmol/L	≤3.0%	≤3.0%	0.2-1.5


*Benchmarking its performance against a gold standard wet bench analyzer



Accuracy Validation

The accuracy of Ansentia Electrolyte Analyzer has been validated by linear regression analysis showing a strong positive, linear relationship between two analyzers ($R^2 > 0.94$ for all parameters), demonstrating the accuracy of our analyzer for K⁺, Na⁺, Cl⁻, iCa²⁺ and iMg²⁺ measurements.

Correlation with Benchtop Analyzer



Test Parameters and Clinical Significance

Potassium (K⁺): Minor fluctuations in extracellular K⁺ concentration can lead to significant alterations in the trans membrane potential gradient, consequently impacting the function of neuromuscular and cardiac tissues.

Sodium (Na⁺): Na⁺ is the principal determinant of water distribution between the intracellular and extracellular compartments, hence Na⁺ is essential for the maintenance of blood volume and thereby blood pressure.

Chloride (Cl⁻): Cl⁻ is essential for the maintenance of normal plasma osmolarity.

Ionized Calcium (iCa²⁺): The maintenance of iCa²⁺ is critical for the structural integrity of bones as well as hemostasis, cardiac and skeletal muscle cell contraction, neuromuscular transmission and the proper function of many hormones.

Ionized magnesium (iMg²⁺): iMg²⁺ is associated with the stabilization of intracellular potassium, which ensures normal functionality of the myocardium, nerves and muscles.

