

# Children's manual wheelchair

## Suitable for:

Individuals with osteogenesis imperfecta (OI), ALS (Amyotrophic Lateral Sclerosis), and penguin-dwelling individuals.

## Function

It can be adjusted at different stages as the child grows. The seat depth, backrest height, and footrest height are adjustable. Accessories such as anti-tipping wheels, push handles, and armrests can be installed or removed, truly achieving customized service.

## Reference data

Material	Frame	Magnesium Alloy
	Flange	Ultralight Carbon Fiber
	Pedals:	Ultralight Carbon Fiber
	Seat	Ultralight Carbon Fiber
	Fork	Aluminum Alloy
	Backrest	Ultralight Carbon Fiber
	Handwheel Rim:	Aluminum Alloy
	Saddle	Sponge
	Spokes	Steel
	Wheels	Aluminum Alloy
	Anti-roll Wheel	Aluminum Alloy
	Rear Push Handle	Aluminum Alloy
	wheel	front wheel
back wheel		24 inches
rear wheel angle		0° , 3° optional
定制	seat width	26 28 30 32 34
	seat depth	27 29 31 33 35 37
	leg length	27 29 31 33 35 37 39 41
	High backrest	30 32 34

	height of seat	Before 50 after42
Add	handrail	aluminum alloy carbon fiber
apolegamy	three-blade wheel	Ultra-lightweight carbon fiber (pair)
	Thickened spokes	steel wire
other	fold	Removable panel
	dismantle	The Great Treatise is Dismantled
	heart	adjustable
	Backrest angle	adjustable
	brake	push brake
	forearm	unadjustable
	seat board	adjustable
	guard board	adjustable
	pedal wheel	adjustable
	front wheel	Adjustable

