

Manual wheel chair

一、Standard

名称	参数
Usage (cm)	65
Folded (cm)	23
Seat Width (cm)	46
Rear Wheel Diameter (cm)	58
Front Wheel Diameter (cm)	20
Seat Cushion Height (cm)	48
Seat Cushion Depth (cm)	40
Overall Length (cm)	104
Overall Height (cm)	87
Backrest Height (cm)	41
Maximum Load (kg)	100
Net Weight (kg)	16
Carton (cm) Unit/Carton	92*22*89

二、Wheelchair performance and structural features:

The product conforms to GB/T13800-2009, the national standard for manual wheelchairs. Its structural features are as follows:

Frame: Constructed from 22mm diameter, 1.2mm thick steel tubing, foldable, with fixed armrests and footrests, electroplated surface, simple, safe, and sturdy structure; Front wheel: High-quality 20cm (8-inch) diameter PVC wheel with electroplated steel fork; flexible and

wear-resistant; Rear wheel: 58cm diameter solid rubber tire with steel wire wheel, electroplated steel handwheel, airless, wear-resistant; Seat and back cushion: Leather fabric, reinforced with canvas and metal strips, and foam for added comfort and easy cleaning;

Brakes: Steel parking assembly, safe and reliable; Armrests: Fixed armrests with long PVC armrest surfaces and stainless steel guards; Footrest: Fixed footrest with high-strength aluminum alloy pedals and PVC leg straps; adjustable and tiltable footrest height.

