


# AST-2876

## Stereomicroscope



目镜 Eyepiece	标准配置 Standard Configuration	
	工作距离 120mm Working Distance	
	放大倍率 Magnification	视场范围(MM) Field of View
WF10X/22	7X	31.4
	60X	3.7
WF20X/10	14X	14.3
	120X	2
放大倍数 Magnification	标准配置: 7X-60X 连续变倍, 选配目镜和辅助物镜可扩展至 3.5X-200X Standard: 7X-60X stereo zoom, 3.5X-200X with optional lens	

<p>目镜 Eyepiece</p>	<p>标准配置: WF10X/22MM 大视野、广角、高眼点, 为佩戴眼睛的观察者提供方便, 可选 WF15X/15MM, WF20X/10MM, WF25X/9MM 目镜</p> <p>Standard WF10X/22MM high eye-point wide field eyepiece, convenience for eye-wearing observers. WF15X/15MM, WF20X/10MM, WF25X/9MM eyepiece is optional.</p>
<p>观察头 Microscope Head</p>	<p>双目/三目观察头, 45°倾斜, 360°可旋转。瞳距调节范围: 54-76mm, 双边视度可调</p> <p>Binocular/Trinocular head, 45°inclined, 360°rotatable. Interpupillary distanc: 54-76mm. Diopter adjustable on both tubes.</p>
<p>物镜 Objective</p>	<p>0.7X-6X 连续变倍, 物镜变倍比 8.6:1</p> <p>0.7X-6X stereo zoom. Zoom Ratio 8.6:1</p>
<p>辅助物镜 Auxiliary Lens</p>	<p>可选配辅助物镜: 0.5X, 工作距离可达到 172mm</p> <p>Can choose 0.5X Auxiliary Lens, the working distance is 172mm</p>
<p>调焦机构 Focusing Adjustment</p>	<p>A1 调焦托架, 手轮松紧可调, 升降范围 50mm</p> <p>A1 focus arm with tension adjustment, rang: 50mm</p>



底座 <b>Microscope Stand</b>	B1 小平板底座，配直径 95MM 黑白工作板，可选配其他底座  B1 stand with black/white plate in diameter 95mm.Other stands are variable.
照明 <b>Illumination</b>	LED 环形灯/光纤照明  LED ring lamp/optical fiber source
三目接口 <b>Camera Adapter</b>	配合三目杆子和 CTV,可接相机、CCD  0.5X,0.35X,1X for option