

# Nano-steam assisted bathing machine

## I、Product presentation

1. Only one person is needed to provide the in-home bathing service;  
No need for water hookups or moving the patient, the service is completed directly in bed, reducing risks;
2. Bathing can be done while clothed, protecting privacy;
3. No need for bath soap, providing deep cleaning and reducing the workload for service personnel;
4. 9 adjustable temperature settings, suitable for bathing all year round;
5. Bathing is done in a sealed environment using a steam canopy on the bed, preventing chills and colds;
6. Steam bathing provides sterilization and deodorization, promotes blood circulation, and strengthens the body.

## II、technical parameter

Input voltage: 220-240V/50HZ	Internal water tank capacity: 1.8 liters
Nano mist output: 40L/min	Power cord length: 1 meter
Power cord length: 1.5 meters	Product weight: 9.8 kg
Packaging box dimensions: 32.5 x 27.5 x 42.5 cm	Rated power: 3000W

