

Portable shower

I、Product presentation

1. The machine operates in a handheld spray shower mode; pulling the trigger on the spray gun starts the spray shower, and releasing the trigger stops it.
2. The portable shower system design eliminates the need to move the assisted person to the bathroom, preventing accidental falls and avoiding unnecessary injuries.
3. The system allows for easy bathing whether the assisted person is sitting or lying down, helping the elderly achieve safe, refreshing, and dignified bathing.
4. The temperature is adjustable: turning clockwise increases heating power, resulting in faster heating and higher temperature; turning counterclockwise decreases power, resulting in slower heating and lower temperature; stopping the knob maintains a constant temperature.
5. Nozzle: Water is sprayed from the nozzle in a nano-atomized form.
6. Environmentally friendly and water-saving design, using approximately 3 liters per use, resulting in low bathing costs.
7. Operating voltage: 220V/single phase; Machine weight: 3KG
8. Total power: 1.5KW; Water consumption: Approximately 3 liters/use
9. Atomization area: Manually adjusted by the handheld atomizer.
10. Adjustable outlet atomization temperature range: Less than 50 degrees Celsius.
11. Machine weight: 3KG
12. Dimensions: 349*141*208mm.

