

Smart shower

Product presentation

1. The machine has three functional modes: body wiping, hair washing, and showering. Each mode has both a standard mode and an enhanced mode.
2. The movable bathing system design eliminates the need to transfer the assisted person to the bathroom, preventing accidental falls and avoiding unnecessary injuries.
3. The assisted person can easily have their body wiped and hair washed whether sitting or lying down. It includes a massage function to help the elderly achieve safe, refreshing, and dignified bathing.
4. To activate the body wiping/hair washing/showering mode, press the clean water outlet control button, and clean water will spray out. Releasing the button will stop the water flow. When all modes (body wiping, hair washing, showering) are off, clean water will not spray out regardless of whether the control button is pressed or released.
5. Built-in high-density 5000mAh lithium-ion power battery allows for operation without being plugged in.
6. The product is equipped with a sponge showerhead and a silicone showerhead. The silicone showerhead is used for washing hair, and the sponge showerhead is used for body wiping and showering.
7. The device has built-in Bluetooth functionality for data interconnection with mobile phones.
8. It uses an integrated washing and suction negative pressure deep cleaning technology. A brushless high-speed fan suction system simultaneously cleans and removes wastewater during bathing, preventing water leakage and cross-contamination, ensuring comfort and hygiene. The powerful suction function, combined with a special bath liquid, deeply cleans and removes dirt from hair follicles, resulting in a fresher feeling after bathing, eliminating body odor, and nourishing hair and skin.

9. The intelligent new assisted bathing device uses recyclable contact points for safe bathing without electric leakage.
10. Environmentally friendly and water-saving design: 3L of water for one hair wash, 6L of water for one body wipe.
11. Ergonomic design, simple and convenient operation, reducing the workload of caregivers and improving service efficiency and quality.
12. Uses a silent brushless high-speed fan and labyrinth noise reduction technology for enhanced comfort for the elderly.
13. The showerhead has a water flow massage function. It allows for direct hand movement of the showerhead for massage and hair washing without wetting the hands, making it convenient and gentle on the hands.
14. The sponge showerhead and silicone showerhead comply with national environmental protection standards;
15. The wastewater tank is built-in, and the clean water tank is external;
16. Equipped with a 4.3-inch display screen to instantly display the machine's working status;
17. The machine has a built-in water heating system, with a water temperature setting range of 22 – 55 °C . Heating stops when the set temperature is reached;
18. The clean water tank has a double-layer insulation function. After power off, the hot water can maintain a temperature of 35 – 45 °C for approximately 4 hours;

