

## KFT-20V Karl Fischer Coulometer



### FEATURES

- LCD display screen.
- Generator cell with diaphragm.
- Selectable units including  $\mu$ g,mg,%,mg/L,g/L,etc.
- Settable parameters,including measurement unit,polarization current,generator current,endpoint potential,stop criterion anetc.
- Data can be stored and easily transferred to printer via RS-232 communication interface.
- Reset feature automatically resumes allsettings back to factor default options.

| Model                |                               | KFT-20V                       |
|----------------------|-------------------------------|-------------------------------|
| Coulometric Ttration | Range                         | 10ug~20mg                     |
|                      | Polarization Current Accuracy | ±0.2μA                        |
|                      | Working Current               | 10mA, 20 mA, 50 mA, 100 mA    |
|                      | Working Current Accuracy      | ±0.5%                         |
|                      | Accuracy                      | ± (5%+3) μg                   |
|                      | Repeatability                 | ≤3%                           |
| General              | Power Supply                  | AC Adapter, 200-240 VAC input |
|                      | Dimensions                    | 300×235×100 mm                |
|                      | Weight                        | 2 kg                          |