

## Tablet Friability Tester



TFT-I



TFT-II/TFT-III



TFT-IV

### Advantage

1. Automation: auto-test, auto-diagnose, auto-alarm.
2. Intelligently control the rotating speed of cylinder and revolutions.
3. One drum runs synchronously. Auto stop at timing point, (only for TFT-IV)
4. Single path and double cylinder, synchronous operation and automatic stop.

5. The revolution of the drum can be preset at any time, either preset or real data is displayed on the LCD.

6. Microprocessor unit is used to control all operations, high accuracy is achieved in testing rotary speed and circle number.

### Specification

Model	TFT-I	TFT-II	TFT-III	TFT-IV
Cylinder Number	2			1
Cylinder Radius	286			
Cylinder Depth	39mm			
Tablet Fall Height	156mm			
Rotating Speed	20~100rpm	25rpm		
Rotation Count	5-900 turns	100 turns	5-900 turns	10-900 turns
Consumption	40W	20W		50W
Power Supply	AC110/220V±10%, 50/60HZ			
Dimensions(WxDxH)mm	370x300x340			300x200x340
Shipping Size(WxDxH)(mm)	460x420x550			430x310x450
N.W./G.W. (kg)	12/20			8/12