

# 20L Automatic Electric Heating Dispenser

**Product Series:** Water Distillers (TT, DZ Series)

**Model:** TT-98-III/DZ

**Internal Volume:** 8.7L / 13.6L /  
25.5L

**Functional Features:** Low  
Water Alarm, Low Water Power-Off, Low  
Water Recovery

## Product Overview:

Distillation of distilled water

## Features

Automatic water level maintenance system.

All-stainless steel design.

Easy operation, stable water output.

DZ type level control, low water alarm, power-off in case of low water.



**Standard: YY/T 0280-1995**

**Safety: Low water alarm, power-off in case of low water.**

## Technical parameters:

Classification	TT-98-III			DZ		
Specifications						
Model	HS.Z11.5L	HS.Z11.10L	HS.Z11.20L	DZ-5LIII	DZ-10LIII	DZ-20LIII
Type	Water vapor condensation	Water vapor condensation	Water vapor condensation	Water vapor condensation	Water vapor condensation	Water vapor condensation
Performance	5	10	20	5	10	20
Condensate Capacity (L/h)						
Components	Stainless steel	Stainless steel	Stainless steel	Stainless steel	Stainless steel	Stainless steel
Evaporator	Stainless steel heating element	Stainless steel heating element	Stainless steel heating element	Stainless steel heating element	Stainless steel heating element	Stainless steel heating element
Heater	4.5kw	7.5kw	15kw	4.5kw	7.5kw	15kw
Rated Power						
Controller	×	×	×	✓	✓	✓
Water Shortage						



Alarm						
Safety Device	φ 305*755	φ 335*855	φ 405*1005	φ 305*755	φ 335*855	φ 405*1005
Dimensions (Diameter * Height mm)	380.35/835	380.36/925	410.41/1135	380.35/835	380.36/925	410.41/1135
Packaging Dimensions (Width * Depth * Height mm)	AC220V/20A	AC380V/11A	AC380V/23A	AC220V/20A	AC380V/11A	AC380V/23A
Power Supply (50/60Hz)	4.7/6.7	6.3/7.8	9.5/11	5/7	7/9	10/12
Rated Current	Stainless steel	Stainless steel	Stainless steel	Stainless steel	Stainless steel	Stainless steel
Net Weight / Gross Weight (kg)	×	×	×	✓	✓	✓

### Automatic water level maintenance system

**Thermostat reset switch (panel indicator: Low water alarm, add water and then turn/press the black switch to reset)**

Automatic water level  
maintenance system



Temperature limiter reset  
switch



Temperature limiter probe

