


Benchtop Low Speed Centrifuge




24× 10ml(4000r/min) RFC:2200×g



6×50ml(4000r/min) RFC:2200×g



12× 10/15ml(4000r/min) RFC:2200×g



18× 10ml(4000r/min) RFC:2200×g

Specification

Model	AST-4LE
Max. Speed	4000r/min
Max. RCF	2200×g
Max.Capacity	6×50ml
Angle Rotor	24×10ml(4000r/min)RFC:2200×g
Angle Rotor	6×50ml(4000r/min)RFC:2200×g
Angle Rotor	12×10/15ml(4000r/min)RFC:2200×g
Angle Rotor	18× 10ml(4000r/min)RFC:2200×g
Time Setting	0~99min
Motor	Brushless DC motor, microcomputer control
Power Supply	AC 220V 50Hz 3A
External Size (WxDxH) (mm)	430x340x250
Shipping Size (WxDxH) (mm)	510x420x330
N.W./G.W. (kg)	20/23