

AST-Gel Card Centrifuge



12 cards swing rotor



24 cards swing rotor

Description

The Gel Card Centrifuge, also named blood group serological test multipurpose centrifuge, it is widely used in the test of blood group serology, bloodtype, microcolumn agglutination, immunodetection and etc.

Advantage

- Brushless motor to make sure free maintenance, no powder pollution, microcomputer control and digital display the time, speed and RCF to let operation easily.
- Special centrifugation technical to different reagent according to the requirement.
- Regular centrifugation program, customers can press the button to choose the required one.
- 10-speed acceleration and deceleration control, booting is the last used program.
- The mute electronic door lock is easy to use and has an emergency door opening function.
- Door lock, overspeed and imbalance automatic protection , Self-diagnostic error.

Specifications

Model	AST-GC12	AST-GC24
Max Speed	4000rpm	3000rpm
Max RCF	1480×g	1192×g
Max Capacity	12 cards	24 cards
Speed Accuracy	±10r/min	±10r/min
Timer	1s~99min59s	1s~99min59s
Noise	≤55dBA	≤55dBA
Power supply	AC110/220V,50/60HZ	AC110/220V,50/60HZ
External Size (W×D×H)	450×320×225mm	500×360×250mm
Shipping Size (W×D×H)	565×410×340mm	615×450×390mm
N.W./G.W.	15kg/17kg	19kg/21kg

