

## Tissue Water Bath



### Advantage

1. Glass tub with adjustable LED illumination with ultra-bright maximum intensity.
2. Water bath can be set from room temperature up to 90°C.
3. Special surface material for durability and corrosion resistance.
4. Provides stable temperature control and buoyancy support.

## Specification

| <b>Model</b>                         | <b>AST-I</b>                                 |
|--------------------------------------|--|
| <b>Temperature range</b>             | <b>RT~90°C</b>                               |
| <b>Temperature control precision</b> | <b>±1°C</b>                                  |
| <b>Dimensions of water bath dish</b> | <b>240×180×50 mm</b>                         |
| <b>Power Supply</b>                  | <b>AC 100V ~ 120V/AC 220V ~ 240V;50/60Hz</b> |
| <b>Consumption</b>                   | <b>200W</b>                                  |
| <b>Dimensions(W×D×H)(mm)</b>         | <b>350×340×120</b>                           |
| <b>Shipping size(W×D×H)(mm)</b>      | <b>430x440x210</b>                           |
| <b>N.W./G.W. (kg)</b>                | <b>6/8</b>                                   |