

Tissue Water Bath



Advantage

1. Large LCD screen with clear display.
2. Unique control microprocessor, maximally meet the needs of different customers.
3. Fast heating efficiency, long life.
4. Energy saving, reliable quality.

Specification

Model	AST-TWB
Temp.setting	RT~90°C (constant adjusting)
Temp.accuracy	±2°C
Ambient temperature	0~40°C
Power Supply	220V/50Hz
Shipping size(W×D×H)(mm)	450x420x200
N.W./G.W. (kg)	7/11