



AST-S7045B



AST-S7045T

## Specifications

Model	AST-S7045B	AST-S7045T
Magnification	Zoom objective 6.7X-45X	
Eyeiece	WF10X/22mm high-eyepoint, wide-field; Convenient for observers wearing glasses *Optional Eyepiece15X/16mm , 20X/12mm , 25X/9mm , 30X/8mm	
Main Body	Binocular head , 45° inclined, 360° rotatable, can be locked in any position desired.  Diopter adjustment:±5 diopters	Trinocular head, 45° inclined, 360° rotatable, can be locked in any position desired. Diopter adjustment:±5 diopters
	Interpupillary Distance adjustment:54-75mm	

Objective	0.67X-4.5X stereo zoom; Zoom Ratio:6.7:1 *Aux Lens: 0.5X , 2X
Focusing Mount	Focusing track mount with vertical working distance 104mm
Stand	Track Stand Vertical Height: 300mm; Mounting Size: $\phi$ 76mm; Focusing Distance: 106mm; With $\phi$ 95mm B/W Plastic & Perspex Plate and Paired Clips; Base Size: 205x275x40mm
Illumination	Top and bottom LED Illumination with an independent control; Input Power : 110v~240v
C-Mount	MSC 1X C-Mount; xMSC 0.5X C-Mount;
Shipping Size(W×D×H)(mm)	510x360x300
G.W. (kg)	9

NOTE: MSC-S7045B 3.35X-270X with optional eyepieces and objectives;

