



SZM-45B



SZM-45T



0.3倍接口 0.5倍接口 1倍接口 0.3倍辅助物镜 0.6倍辅助物镜 0.7倍辅助物镜 0.75倍辅助物镜 1.5倍辅助物镜 2倍辅助物镜

放大倍数	标准配置: 7X-45X,连续变倍 Zoom objective
Magnification	7X-45X
目镜Eyepiece	WF10X/20mm,大视野、广角、高眼点、为佩戴眼镜的观察者提供方便 WF10X/20mm high-eyepoint,wide-field;Convenientfor observers wearing glasses *目镜Optional Eyepiece15X/15mm,20×/10mm,25X/9mm,30X/8mm
观察头	铰链式双目/铰链式三目, 45°倾斜、360°旋转, 瞳距54-76mm;双边视度调节(±5)
Main Body	Binocular head/Trinocular head,45°inclined,360°rotatable,can be locked in any position desired.Interpupillary Distance adjustment54-76mm;Dioptr adjustment:±5 diopters
物镜Objective	0.7X-4.5X连续变倍、物镜变信比6.4:10.7X-4.5X stereo zoom;Zoom Ratio:6.4:1 *辅助物镜Aux Lens:0.5X,0.7X,0.75X,1.5X,2X
调焦托架	调焦托架手轮松紧可调、升降范围50mm
Focusing Mount	Focusing mount with vertical working distance 50mm

底座 Stand	立杆卤素灯底座 立柱高：240mm;直径：φ32mm;配有压片和φ95mm黑白板和玻璃板；底座： 200×255x60mm Vertical Post:240mm;Diameter:φ32mm;With φ95mm B/W Plastic &Perspex Plate and Paired Clips;Base Size:200×255×60mm
照明 Illumination	配上下光源12v15w卤素灯； *LED光源， *7w荧光灯或者环形灯 Transmitted and reflected 12v/15w halogen;*7w fluorescent,*LEDillumination
三目接口 C-Mount	1倍接口， *0.5倍接口， *0.3倍接口 SZM 1XC-Mount;*SZM 0.5XC-Mount;*SZM 0.3XC-Mount
注：选配目镜和辅助物镜可扩展至3.5X-270X,连续变倍 Information:3.5X-270Xwith optional eyepieces and objectives *为选配Optional configuration	

