



济南安生泰医疗器械有限公司
JINAN ANSENTA MEDICAL DEVICES CO., LTD.

Add:Room 1208,Century Fortune
Center Building A,Sanqing,Jinan
City
Tel: + 86 17661957135
Email:info@ansenta.com

NG-300 High-Purity Nitrogen Generator

- Gas is produced simply by turning on the switch; it is safe and reliable.
- The output flow rate can be automatically adjusted according to the usage.
- Equipped with overpressure protection and anti-overflow device
- LED digital display
- Impact-resistant and aging-resistant filter
- High gas purity with no purity reduction
- All gas lines are made of high-quality stainless steel.
- Long service life and easy maintenance
- Can be used 24 hours a day



It can be used with gas chromatographs of any model and manufacturer, both domestically and internationally.

Add alkali once; only distilled water needs to be replenished for routine maintenance.

Nitrogen generator technical parameters:

Product Model	NG-300
Output flow	0~300ml/min
Work pressure	0~0.4Mpa
Nitrogen purity	99.999% oxygen content<3PPM,Water-containing dew point-56°C
Operating voltage	~220V±10% 50Hz
Ambient temperature	5°C~40°C
Rated power	150W
External dimensions	400×200×300 (mm)
Instrument weight	8kg