

## Precision Washing Elisa Plate Washer



### Advantage

- 1. Large LCD with menu design for easy use;
- 2. Whole plate cleaning or single strip cleaning;
- 3. Waste liquid overflow alarm function;

4. One-point, two-point and whole bottom washing;
5. Flat bottom, V-bottom or U-bottom strip cleaning;
6. Large memory for storing up to 48 user-programmed cleaning programs.

## Specification

Model	ELMW-310	ELMW-320
Input	8 key membrane keyboard	
Display	High brightness LCD display screen	
Manifold	8 pins and 12 pins	
Wash Mode	Strip mode and plate mode	
Wash strips	1~12 adjustable	
Applicable well	Flat, V-bottom or U-bottom	
Washing times	0~99 times adjustable	
Dipping time	0~3600s adjustable	
Shaking time	0~600s adjustable	
Residual volume	$\leq 1\mu\text{L}/\text{well}$	

Liquid volume	50~3000ul/well	
Sipping time	0.1~9.9s	
Storage Capacity	More than 100 user defined wash board procedures	
Power Supply	220V,50Hz	
Washing fluid channel	1	1(Optional: 2-4 Channels)
Dimensions(W×D×H)(mm)	448x382x163	
Shipping size(W×D×H)(mm)	620x550x490	
N.W./G.W. (kg)	10/16	