

Tissue Water Bath



Advantage

1. Large LCD screen with clear display.
2. Unique control microprocessor, maximally meet the needs of different customers.
3. Fast heating efficiency, long life.
4. Energy saving, reliable quality.

Specification

Model	AST-TWB90
Temperature range	RT~90℃
Temperature control precision	±1℃
Dimensions of water bath dish	260×220×65mm
Material of internal bath	Aluminum
Power Supply	AC220V 50Hz ; AC110V 60Hz;50/60Hz
Consumption	400VA
Fuse	F6A L 250AC
Dimensions(W×D×H)(mm)	340x400x140
Shipping size(W×D×H)(mm)	550x490x260
N.W./G.W. (kg)	6/8

Tissue Water Bath



Advantage

1. Large LCD screen with clear display.
2. Unique control microprocessor, maximally meet the needs of different customers.
3. Fast heating efficiency, long life. 4. Energy saving, reliable quality.

Specification

Model	AST-TWB85
Temperature range	RT~85℃
Temperature control precision	±1℃
Dimensions of water bath dish	350×220×45mm
Power Supply	AC 100V~120V/AC 220V~240V;50/60Hz
Consumption	400W
Dimensions(W×D×H)(mm)	510×420×150
Shipping size(W×D×H)(mm)	620x520x270
N.W./G.W. (kg)	9/15