

## Ultrasonic Height and Weight Analyzer



### Advantage

- Ultrasonic height.
- Foldable, with wheels at the bottom.
- With the RS232 port, data can be transferred to the computer.
- With a printer function.

## Specification

Model	UHWA-100
Height Measurement	Ultrasonic intelligent induction measurement
Height Measurement range	50-200cm,d=0.5cm
Weighing Range	3-200kg,d=0.1kg
Measurement Items	Weight,Height,BMI
Voice Reading Function	voice can be on or off
Interface	With RS232 Interface
Printing Function	comes with a thermal printer
Overweight Display	Err; Loaded weight display: 0.0kg
LCD large Screen,	size is 67×77mm
Switching Mode	Gravity sensing power on/power-on button ; auto power off/power-off button
Power Supply	Plug-in use ( equipped with a DC 6V2A adapter)
Packing Size	120×45×27.5cm
Net Weight	16.7kg
Gross Weight	23kg

