

Vertical Pressure Steam Sterillzer

Digital Display Self-control



Advantage

- Hand wheel translation and quick opening structure.
- High quality stainless steel material.
- Microcomputer automatic control, arbitrary setting. sterilization parameters.
- Safety interlocking device.

- Digital display and touch buttons.
- Standard test interface.
- Self-expanding seal.
- With automatic protection function: over temperature protection; over pressure self discharge protection, low water level protection; anti dry burning.
- With stainless steel bucket or basket.
- Buzzer reminder after sterilization, automatic stop.
- Automatic discharge of cold air, automatic steam exhaust after sterilization.

Specification

Model	ST-VB35	ST-VB50	ST-VB75	ST-VB100
Capacity (L)	35	50	75	100
Consumption (kW)	3.5	3.5	4.5	4.5
Power Supply	AC 220/110V, 50/60Hz			
Design Pressure	0.25MPa			
Design Temperature	139°C			
Rated Working Pressure	0.22MPa			
Rated Working Temperature	134°C			
Sterilization Temperature Setting Range	116°C~134°C			
Sterilization Time Setting Range	4~120min			
Chamber Capacity (Ø×H) (mm)	300×450	350×515	400×585	450×615
Bucket Size (Ø×H) (mm)	275×395	325×235×2	375×265×2	425×275×2
Basket Size (Ø×H) (mm)	270×380	320×220×2	370×250×2	420×265×2
External Size (W×D×H) (mm)	476×473×990	544×536×1100	586×576×1210	620×610×1275
Shipping Size (W×D×H) (mm)	630×580×1130	630×580×1230	680×630×1290	650×690×1410

Model	ST-VB35	ST-VB50	ST-VB75	ST-VB100
Gross Weight (kg)	79	106	120	133
Net Weight (kg)	58	80	90	106