

HH 系列恒温水浴锅（新）



Product Overview:

The inner liner of the constant temperature water bath is made of high-quality 304 stainless steel, while the outer shell is formed by stamping, bending, and welding cold-rolled steel. The outer shell surface is coated with a high-strength electrostatic powder coating, resulting in a smooth and durable paint film. It can be used for drying, concentration, distillation, impregnation of chemical reagents, pharmaceuticals, and biological products. It can also be used for constant temperature heating in a water bath and other temperature experiments. It is an essential tool for biological, genetic, virological, aquatic, environmental, pharmaceutical, health, biochemical, analytical, and educational research laboratories.

Product Features

1. Utilizes LED digital displays with an elongated instrument design, standard timer function, clear and intuitive operation, and PID temperature control technology for rapid heating and accurate temperature control.



济南安生泰医疗器械有限公司
JINAN ANSENTA MEDICAL DEVICES

Add:Room 1208,Century Fortune
Center Building A,Sanqing,Jinan
City
Tel: + 86 17661957135
Email:info@ansenta.com
www.ansenta.com

2. Cold-rolled steel outer shell with electrostatic powder coating, sturdy, durable, and corrosion-resistant.

3. 304 stainless steel inner liner, lid, and shelf, corrosion-resistant, bright, and easy to clean.

4. Ring-shaped top cover, adaptable to various container sizes, improving heat preservation performance.

5. Equipped with a magnetic stirring function to improve water bath temperature uniformity (AST-HH-8 only).

6. Built-in rubber drain pipe for convenient drainage (AST-HH-4/6/8 only).

Technical parameters

Model Number	AST-HH-1	AST-HH-2	AST-HH-4	AST-HH-6	AST-HH-8
Power Supply Voltage	AC220V 50HZ				
Power (W)	350	600	1000	1500	1800
Number of Holes	1	2	4	6	8
Display Method	LED Digital Tube				
Temperature Control Range	RT+5-100°C				
Temperature Control Accuracy	±1°C			±0.5°C	
Temperature Display Accuracy	0.1°C				
Timing Range	0-999min				
Working Dimensions (L)	160×160×105	305×160×105	305×305×105	470×305×105	610×305×105



济南安生泰医疗器械有限公司
JINAN ANSENTA MEDICAL DEVICES

Add:Room 1208,Century Fortune
Center Building A,Sanqing,Jinan
City
Tel: + 86 17661957135
Email:info@ansenta.com
www.ansenta.com

×W×H) mm					
Inner Tank Dimensions (L ×W×H) mm	160×160×140	305×160×140	305×305×140	470×305×140	610×305×140
Net Weight (kg)	2.6	3.8	5.3	7.5	9.8
Product Dimensions (L ×W×H) mm	200×200×200	345×200×200	345×345×215	510×345×240	650×345×260

Note: LC-HH-8 has a stirring function for water bath media.

No.	Description	Quantity
1	Main unit (including sample rack and cover plate)	1 unit
2	Ring-shaped pot cover	Determined by number of holes
3	Power cord	1 pc
4	Fuse	1 pc
5	Magnetic stir bar (AST-HH-8)	1 pc
6	Instruction manual	1 copy
7	Certificate of conformity / Warranty card	1 copy

Packing list