



TDZ4-WS Benchtop Low Speed Centrifuge

TDZ4-WS benchtop low speed centrifuge is equipped with four angle rotors: 18x10ml, 12x20ml, 4x50ml 24x10ml and one swing rotor 8x10ml . It is widely used in most of the hospitals and universities by its high performance.



Product advantages

- ◆Microprocessor control, DC brush-less motor.
- ◆Touch panel and LCD display.
- ◆Special damping structure to reduce vibration.
- ◆Countdown time displayed by seconds if the centrifugation time is less than 1min.



Technical parameter

| | |
|---------------------------|--------------------|
| Order NO | 190000 |
| Max Speed | 4000r/min |
| Max RCF | 2200xg |
| Max Capacity | 12×20ml |
| Speed Accuracy | ±30r/min |
| Time Setting Range | 1min~99min |
| Noise | ≤ 62dB (A) |
| Power Supply | AC 220±22V 50/60Hz |
| Power Consumption | 100W |
| Dimension (WxDxH) | 330×390×250(mm) |
| Packing Dimension (WxDxH) | 395×475×290 (mm) |
| Net Weight | 16kg |



Rotor parameters



Rotor Type: Angle rotor
Max Speed: 4000r/min
Max RCF: 2200xg
Max Capacity: 12×20ml
Order NO: 0301808500
Adapter: 5ml, 7ml Blood collection tube



Rotor Type: Angle rotor
Max Speed: 4000r/min
Max RCF: 2200xg
Max Capacity: 18×10ml
Order NO: 0301808400
Adapter: 5ml, 7ml Blood collection tube



Rotor Type: Angle rotor
Max Speed: 4000r/min
Max RCF: 2200xg
Max Capacity: 24×10ml
Order NO: 0301808300
Adapter: 5ml, 7ml Blood collection tube



Rotor Type: Angle rotor
Max Speed: 4000r/min
Max RCF: 2200xg
Max Capacity: 4×50ml
Order NO: 0301808600